

Thank you for your purchase of Beam Professional Teeth Whitening products. All of our whitening solutions are manufactured in a FDA registered facility in the USA.

## We Provide Effective And Affordable Teeth Whitening Solutions Using The Highest Quality Ingredients Available.

### STEP 1

Rinse the LED Smart Tray with warm water.

### STEP 2

Apply whitening solution to your teeth. Rotate the pen at the base to dispense the gel onto the brush. Make a big wide smile and apply the gel generously to the teeth using the brush. Be sure to apply the whitening gel to the front of the teeth only.

### STEP 3

Connect the LED Smart Tray to any smart device or power bank of your choice and turn it on so that the LED light is constant.

### STEP 4

Put the Smart Tray into your mouth and bite down firmly. Keep the tray in your mouth for 20-30 minutes.

### STEP 5

Remove the tray from your mouth and rinse your mouth and tray with lukewarm water and allow the tray to dry before storing it.

**There are up to 15 treatments included with one pen.**

We recommend when using the system for the first time to use it daily until desired results are achieved. Then reduce to 1-2 times per week to maintain.

To get more whitening solution pens please go to

[www.BeamTeethWhitening.com](http://www.BeamTeethWhitening.com)

### DISCLAIMERS:

Keep out of reach of children. Avoid contact with eyes and skin. Consult your dentist before using. Do not use if pregnant, lactating, have gum disease or other dental conditions. Discontinue use if pain on teeth or gums occurs. Use only as instructed. Use at your own risk. Purchasing any of our products frees all distributors, manufacturers and affiliates from any liability associated with the use of this system. Do not use if you are allergic to peroxide products. The whitening solution is a 12% Hydrogen Peroxide whitening solution. 12% HP is equivalent to 36% Carbamide Peroxide. A temporary side effect of tooth sensitivity or mild tissue irritation is possible but usually transient and may last for a day or two. Avoid beverages that contain citrus and hot and cold drinks for the next 24 hours after treatment. Sticking to room temperature drinks is ideal. If you have any adverse reactions terminate the use of this product and consult a healthcare professional.