

CARE diagnostica®

viola®

viola®
Pregnancy Test

In-vitro diagnostic medical device
for self-testing
Contents: instructions for use, foil pack
with test stick and desiccant



•reliable
•simple
•fast



To be applied
from the **FIRST** day
of the missed period.

CARE diagnostica
Produktions- und
Vertriebsges.m.b.H
Roemerstrasse 8
2514 Traiskirchen
Austria
www.carediagnostica.at

viola® Pregnancy Test

CE
0483

**Pregnancy
Test**



Keep out of reach
of children.

For external
use only.

Store at
4-30 °C.



Contents of the package: 1 test

Reliable

viola® is being successfully used in clinics
and surgeries.

Simple

Expose the tip of the test stick to the
urine flow for a few seconds.

Fast

Read results after a few minutes.

Applicable at any time of the day.

REF
970110-E

P0004-GB
Rev03
2021-03



viola®
Pregnancy Test