



viola[®] Super Early Pregnancy Test
from day 8

CARE diagnostica[®]

viola[®]

viola[®] Super Early Pregnancy Test
from day 8




Super Early Pregnancy Test
from day 8

2
TESTS

CARE diagnostica
Produktions- und
Vertriebsges. m.b.H.
Roemerstrasse 8
2514 Traiskirchen
AUSTRIA
www.carediagnostica.at

Keep out of reach
of children

For external use only



CE
0483

Clear results from 8 days after conception.

This test is twice as sensitive as 'early' pregnancy tests,
giving you a reliable result in just 10 minutes.

You can use this test even before your period is due -
from 8 days after conception.

Early morning urine test- from 5 IU hCG/L.

In-vitro diagnostic medical device for self-testing

Contents: Instruction leaflet, 2 foil packs each
containing a test stick and desiccant

P0065-GB
Rev04 2023-07

REF
973005-E