



# ***Reusable Cannulas Quick Reference Catalog***

- ***Air Injection/Irrigation***
- ***Alpha-Chymotrypsin***
- ***Capsule Polishers***
- ***Chamber Maintainers***
- ***Cortex Extraction***
- ***Cyclodialysis***
- ***Cystotomes***
- ***Hydrodissection***
- ***I & A Cannulas***
- ***Irrigating Vectis***
- ***Irrigating Choppers***
- ***Lacrimal***
- ***LASIK & LASEK***
- ***Miscellaneous***
- ***Retrobulbar,  
Peribulbar  
& Subtenon***



---

**ELEVATING THE QUALITY OF CARE IN OPHTHALMOLOGY**  
**800-637-4346**

# CONTENTS

<b>About The Catalog</b> .....	1
<b>Air Injection/Irrigation</b> .....	2
<b>Alpha-Chymotrypsin</b> .....	5
<b>Chamber Maintainers</b> .....	5
<b>Capsule Polishers</b> .....	6
<b>Cortex Extraction</b> .....	8
<b>Cyclodialysis</b> .....	10
<b>Cystotome</b> .....	10
<b>Hydrodissection</b> .....	12
<b>I &amp; A Cannulas</b> .....	14
<b>Irrigating Vectis</b> .....	15
<b>Irrigating Choppers</b> .....	16
<b>Lacrimal</b> .....	17
<b>LASIK &amp; LASEK</b> .....	17
<b>Miscellaneous</b> .....	20
<b>Retrobulbar Peribulbar &amp; Subtenon</b> .....	23
<b>Index</b> .....	24

## About Cannulas

Most Of Our Irrigating And Aspirating Cannulas Have Irrigation On The Tubing Hubs And Aspiration On The Needle Hubs. The Cannulas Have 304 Series Stainless Steel Working Tips And The Hubs Are 300 Series Stainless. Each Cannula Undergoes A Rigorous Qualification Process To Ensure Proper Functionality With No Bio-Compatibility Consequences. All Cannulas Can Be Disinfected Or Autoclaved Using Common Sterilization Protocol.

The Following Symbols Are Illustrated Throughout The Section For Quick Function Identification:



<b>GAUGE CONVERSION CHART</b>				
Gauge	ID* (inch)	ID* (mm)	OD** (inch)	OD** (mm)
19g	.028	.70	.043	1.10
20g	.024	.60	.036	.90
21g	.020	.51	.032	.80
22g	.016	.41	.028	.70
23g	.013	.33	.025	.60
24g	.012	.30	.022	.55
25g	.010	.25	.020	.50
26g	.010	.25	.018	.45
27g	.008	.20	.016	.40
28g	.007	.18	.014	.36
29g	.007	.18	.013	.33
30g	.006	.15	.012	.30
<b>*Inner Diameter</b>		<b>**Outer Diameter</b>		





## ***About the Catalog***

All Photographs, Illustrations, And Text In This Catalog Are The Property Of Rhein® Medical, Inc. They Are Copyrighted, And May Not Be Reproduced Without The Written Permission Of Rhein Medical.

Most Of Our Proprietary Instruments May Have A Video On Our Website: [www.RheinMedical.com](http://www.RheinMedical.com)

## ***General Company Information***


We Employ People At Rhein Medical, Not Answering Machines. When You Call Us During Our Normal Business Hours, You Will Reach A Responsive Customer Service Professional, Not An Electronic Routing System. We Will Answer Your Questions, Or Process Your Orders, With A Minimum Of Delay. Our Most Important Concern Is You. We Operate Customer Service Hours From 8:30 AM To 5:00 PM E.S.T., And Orders And Questions May Be Presented Monday Through Friday During These Hours. We Offer A Toll Free Telephone Number In The U.S.A.: 800-637-4346. Orders Or Questions May Also Be Presented By Facsimile: 727-341-8123, Or By Email: [Info@RheinMedical.com](mailto:Info@RheinMedical.com). International Calls May Be Placed By Calling Our Main Telephone Number: 727-209-2244, Or Anyone May Live Chat With An Operator at [www.RheinMedical.com](http://www.RheinMedical.com).

You May Also Place Your Questions And Orders Through Our International Distribution System. Please Contact Rhein Medical Or Check Our Website For The Distributor In Your Territory Or Country.



[www.RheinMedical.com](http://www.RheinMedical.com)



**A**

**B**

**C**



### ***Air Injection/ Irrigation Cannulas***

- A)** 91-1000 Delron Irrigating Handle, Two Male Fittings.
- B)** 91-1001 Delron Irrigating Handle, Male/Female Fittings.
- C)** 91-1002 Delron Irrigating Handle, Two Male Fittings, 75mm Length.
- D)** 91-1003 Aluminum Irrigating Handle, Male/Female Fitting.
- E)** 91-1100 Silicone Squeeze Bulb
- F)** 91-1101 Bishop-Harmon Adaptor For Squeeze Bulb.
- G)** 91-1105 Bishop-Harmon Cannula, W/Bulb & Adaptor.
- H)** 91-2100 Wilbrandt-Binkhorst, 23g, Straight
- 91-2101 Wilbrandt-Binkhorst, 23g, Left

**D**

**E**

**F**

**G**
**H**






**A**

**B**

**C**







### ***Air Injection/ Irrigation Cannulas***

- A)** 91-2102 Wilbrandt-Binkhorst Cannula, 23g, Right.
- B)** 91-2103 Gimbel Irrigating Cannula, 30g, 28mm.
- C)** 91-2104 Gimbel Irrigating Cannula, 27g, 25mm.
- D)** 91-2105 Gills/Welsh Cannula, 25g Angled 9mm Tip.
- E)** 91-7024 Air Injection Cannula, 25g, 5mm Tip.
- F)** 91-7026 Air Injection Cannula, 27g, 7mm Tip.
- 91-7025 Air Injection Cannula, 27g, 5mm Tip.
- G)** 91-7027 Air Injection Cannula, 30g, 7mm Tip.
- 91-7030 Air Injection Cannula, 30g, 5mm Tip

**D**

**E**

**F**

**G**


**A**

**B**

**C**

**D**

**E**

**F**

**G**

**H**


### ***Air Injection/ Irrigation Cannulas***

- A)** 91-7035 Bishop-Harmon Cannula, 19g, 7mm Tip, 24mm.
- B)** 91-7036 Gills I/A Cannula, 25g Aspirating & 23g Irrigating.
- C)** 91-7037 Gills/Welsh Cannula, 25g, Angled 5mm Tip.  
91-7038 Gills/Welsh Cannula, 25g, Angled 8mm Tip.
- D)** 91-7039 Gills/Welsh Cannula, 25g, Angled 10mm Tip.  
91-7040 Gills/Welsh Cannula, 25g, Angled 7mm Tip.
- E)** 91-7097 Knolle Anterior Chamber Cannula, 23g, Angled 2mm Tip, 20mm.  
91-7098 Knolle Anterior Chamber Cannula, 25g, Angled 2mm Tip, 20mm.  
91-7099 Knolle Anterior Chamber Cannula, 27g, Angled 2mm Tip, 20mm.  
91-7103 Knolle Anterior Chamber Cannula, 30g, Angled 2mm Tip, 20mm.
- F)** 91-7127 Anterior Chamber Cannula, 23g Angled, 12mm.  
91-7127-S Straight Anterior Chamber Cannula, 23g , Length 17mm.

**G)** 91-7128 Anterior Chamber Cannula, 26g Angled, 12mm.

**H)** 91-8046 Gimbel Irrigating Cannula, Blunt Tip Angle to Side End Open 26g.


**A**

**B**

### ***Air Injection/ Irrigation Cannulas***

**A)** 91-8042 Huang Irrigation Cannula.


### ***Alpha-Chymotrypsin***

**B)** 91-7045 Troutman Alpha-Chymotrypsin Cannula, 25g, 28mm Length.

**C)** 91-7047 Roper Alpha-Chymotrypsin Cannula, 19g.

91-7048 Roper Alpha-Chymotrypsin Cannula, 19g, Right.

91-7049 Roper Alpha-Chymotrypsin Cannula, 19g, Left.


**C**

**D**

### ***Chamber Maintainers***



**D)** 91-7100 Heslin Infusion Cannula, 23g, 9mm Sandblasted Tip.

91-7095 Heslin Infusion Cannula, 25g, 9mm Sandblasted Tip.

**E)** 91-8013 McIntyre III Chamber Maintainer, Cannula, 19g Tip.

**F)** 91-8037 Steinert Threaded 19g Chamber Maintainer With Tapered Tip.

**G)** 91-8048 Gills Chamber Maintainer, 21g Tip.





**E**

**F**

**G**

**A**


**B**


### ***Capsule Polishers***

- A)** 91-7101 Brown Capsule Polisher, Cannula, 26g.
- B)** 91-7090 Kratz Capsule Scraper, Cannula, 23g, 1mm Sandblasted Tip.
- C)** 91-7091 Knolle Capsule Scraper, 23g, 2mm Sandblasted Tip
- D)** 91-7091-1 Knolle Capsule Scraper, Cannula, 23g, 2mm Sandblasted Tip.
- E)** 91-7091-2 Lessner Capsule Scraper 27g, 3mm Tip, Sandblasted 60° Top & Bottom.
- F)** 91-7092 Jensen Capsule Scraper, Cannula, 27g, Sandblasted Tip.
- G)** 91-7093 Simcoe Capsule Scraper, Cannula, 23g, 13mm Length.
- H)** 91-7093R Simcoe Reverse Aperture Posterior Scrubber, 23g 13mm length.

**C**

**D**

**E**

**F**


**G**

**H**




**A**

**B**

**C**

**D**

**E**

**F**

**G**

**H**

## **Capsule Polishers**

- A)** 91-7094 Draws Capsule Polisher, Cannula, 25g, Port At 6 O'Clock.
- B)** 91-7096 Kraff Capsule Polisher, Cannula, 25g Irrigation Port/Center.
- C)** 91-7101S Brown Capsule Polisher, 26g W/Sandblasted Tip.
- D)** 91-7105 Graether Collar Button, Cannula, 23g Angled.
- 91-7106 Graether Collar Button, Cannula, 23g, Straight.
- E)** 91-8002 Knolle Capsule Scraper, 30g Cannula, 2mm Sandblasted Tip.
- F)** 91-8003 Jensen Capsule Scraper, 25g Cannula, Sandblasted Tip.
- G)** 91-8004 Jensen Capsule Scraper, 27g Cannula, Semi-Sandblasted Tip.
- 91-8005 Jensen Capsule Scraper, 23g Cannula, Sandblasted Tip.
- H)** 91-8006 Jensen Capsule Scraper, 25g Cannula, Semi-Sandblasted Tip.
- 91-8007 Jensen Capsule Scraper 23g Cannula, Semi-Sandblasted Tip.

**A**

**B**


### **Cortex Extraction**

**A)** 91-7041 Simcoe Cortex Extractor, Cannula, 23g, 0.3mm Port, 13mm Tip.

**B)** 91-7042 Simcoe Cortex Extractor, Cannula, 23g, 0.3mm Port, 21mm Tip.

91-7043 Welsh Flat Olive Cortex Extractor, Canula, 25g, Right.

**C)** 91-7044 Welsh Flat Olive Cortex Extractor, Cannula, 25g, Left.


**D)** 91-7046 Welsh Flat Olive Cortex Extractor, Cannula, 25g, Straight.


91-7058 Simcoe Cortex Extractor, Cannula, 23g, 0.4mm Port, 21mm.

**E)** 91-7061 Binkhorst-McIntyre, 26g, Cannula, 18mm Right.

**C**

**D**

**E**



**A**

**B**

**C**

### **Cortex Extraction**

**A)** 91-7062 Binkhorst/McIntyre, 26g, Cannula, 18mm Left.

**B)** 91-7067 Binkhorst-McIntyre, 26g, Cannula, 18mm Straight.

**C)** 91-7068 Binkhorst Cannula, 22g Right, 18mm Length.

**D)** 91-7069 Binkhorst Cannula, 22g Left, 18mm Length.

**E)** 91-7110 J-Shaped Cannula, 21g.

91-7111 J-Shaped Cannula, 23g.

91-7112 J-Shaped Cannula, 25g.


**F)** 91-7114 J-Shaped Cannula, 26g.

**G)** 91-8045 Anis Cortex Aspirating Cannula, Straight, 0.4mm Sideport, 23g Thin Wall.

**H)** 91-8047 Gimbel Crtex Aspirating Cannula End Open 30g.


**D**

**E**

**F**

**G**

**H**

**A**

**B**

**C**


### **Cyclodialysis**


- A)** 91-7050 Kelman Cyclodialysis Cannula, 19g, Flat Right.
- B)** 91-7051 Kelman Cyclodialysis Cannula, 19g, Flat Left.
- C)** 91-7052 Randolph Cyclodialysis Cannula, 19g, Front Opening.
- 91-7053 Randolph Cyclodialysis Cannula, 21g, Front Opening.

**D**


### **Cystotome**


- D)** 91-7001 Formed Cystotome Cannula, 12mm From Tip, 22g.
- 91-7002 Formed Cystotome Cannula, 12mm From Tip, 25g.
- 91-7003 Formed Cystotome Cannula, 12mm From Tip, 27g.
- E)** 91-7004 Drews Cystotome Cannula, Formed, 30g, 15mm Length.
- F)** 91-7010 Kelman Cystotome Cannula, 22g Sharp, Straight.
- 91-7011 Cystotome Cannula, 24g, Straight.
- G)** 91-7012 Knolle-Kelman Cystotome Cannula, 22g Long Sharp, 32mm Length.

**E**

**F**

**G**


**A**

**B**

**C**

**D**

**E**

**F**



**G**


### **Cystotome**


- A)** 91-7015 Lippman Micro Cystotome Cannula, 22g, 23mm.
- B)** 91-7016 Kratz Cystotome Cannula, 22g, Angled 60 Degrees.
- 91-7018 Kratz/Kelman Cystotome Cannula, 22g, Angled, 45 Degrees.
- C)** 91-7019 Kratz/Kelman Formed Cystotome Cannula, 22g
- D)** 91-7020 Nevyas Double Sharp Cystotome Cannula, 22g, 32mm.
- E)** 91-7021 Nevyas Double Sharp Cannula, 22g, 23mm Length, Straight.
- F)** 91-7022 Kelman Cystotome Cannula, Double Sharp Front & Back, 22g.
- G)** 91-7023 Reverse Cutting Cystotome, 23g, Angled 45 Degrees.

**A**

**B**


**C**

**D**

**E**


**F**

**G**


### **Hydrodissection**

- A)** 91-7129 Hydrodissection Cannula, J-Shaped, Flattened, 30g, Straight.
- B)** 91-7131 Hydrodissection Cannula J-Shaped, Flattened, 23g, Straight.
- C)** 91-7132-A Hydrodissection Flattened Cannula, 23g, Angled.
- D)** 91-7132-S Hydrodissection Flattened Cannula, 23g, Straight.
- E)** 91-7133 Hydrodissection Cannula, J-Shaped, Flattened, 27g, Straight.
- F)** 91-7133-A Hydrodissection Cannula J-Shaped 27g Flattened Angled Shaft.
- G)** 91-7134-S Hydrodissection Cannula, 25g, Straight, Straight.
- 91-7134-A Hydrodissection Flattened Cannula, 25g, Angled.


**A**

**B**

**C**

### **Hydrodissection**

**A)** 91-7135-A Hydrodissection Flattened Cannula, 27g, Angled.

**B)** 91-7135-S Hydrodissection Cannula, 27g, Straight.

**C)** 91-7136-A Hydrodissection Flattened Cannula, 30g, Angled.

91-7136-S Hydrodissection Flattened Cannula, 30g, Straight.


**D)** 91-7137 Hydrodissection Cannula, J-Shaped, Flattened, 25g, Straight.

**E)** 91-7138 Kellan Hydro-Elevation 27g Cannula.

**F)** 91-7141 Khouri Hydrodissection Cannula, 23g, Left.

91-7142 Khouri Hydrodissection Cannula, 23g, Right.

**G)** 91-7206 Utrata 27g Hydrodissection "J" Cannula.




**D**

**E**

**F**

**G**

**A**

**B**

**C**


### ***Hydrodissection***

- A)** 91-8011 McIntyre III 27g Hydrodissection Cannula.  
**B)** 91-8012 McIntyre III 21g Debulking Cannula.  
**C)** 91-8016 Ruminson Hydrodissection Luxation Cannula  
**D)** 91-8032 Chang 27g Hydrodissection Cannula, Short 90 Degree Beveled Tip  
 91-8032-L Chang 27g Hydrodissection Cannula, Long 90 Degree Beveled Tip

**D**


### ***Irrigating & Aspiration***

**E**


- E)** 91-7036 Gills I/A Cannula, 25g Aspirating & 23g Irrigating.  
**F)** 91-7055 Simcoe Double Barreled I/A Cannula, 23g Thin Wall.  
**G)** 91-7057 Simcoe Nucleus Spatula, Cannula, 25g, Flattened Tip.  
**H)** 91-7060 Simcoe Double Barreled I/A Reverse Cannula, 23g Thin Wall.

**F**

**G**

**H**


**A**

**B**


### ***Irrigating & Aspiration***

- A)** 91-7063 Gills I/A, 25g Thin Wall/Aspiration,  
**B)** 91-7064 J-Shaped I/A Cannula, 26g Aspiration & 23g Irrigation.  
**C)** 91-7071 Grandon I/A Cannula, "03", 23g X 23g, Asp. Ports TopLeft  
 91-7072 Grandon I/A Cannula, "04", 23g X 23g, Asp. Ports Top Right.  
**D)** 91-7122 McIntyre II 3mm Irrigating Removal Spoon, Small.  
**E)** 91-7144 Khouri J-Shaped I/A Cannula, 26g Irrigation, 23g Aspiration.

**C**


**D**




### ***Irrigating Vectis***

- F)** 91-7065 Knolle/Pearce Vectis Cannula, 23g, 3 Outer Ports.  
 91-7066 Sheets Irrigation Vectis Cannula, 23g Without Serrations.  
**G)** 91-7066-S Sheets Irrigation Vectis Cannula, 23g With Serrations.

**E**

**F**

**G**





### ***Irrigating Vectis***

- A)** 91-7070 Simcoe Nucleus Loop Left, Cannula, 25g, Two Ports, 4mm.
- B)** 91-7075 Simcoe Nucleus Loop Right, Cannula, 25g, Two Ports, 4mm.
- C)** 91-7076 Lens Loop, Cannula, 3mm Wide, 3 Port.
- D)** 91-7123 McIntyre II 4.5mm Irrigating Removal Spoon, Cannula, Large.
- E)** 91-7143 Khouri 3mm Nucleus Removal Spoon, Cannula.

### ***Irrigating Choppers***

- F)** 91-7145 Gillum Irrigating Phaco Splitter.
- G)** 91-7146 Olson 21g Irrigating Quick Chopper, Right Hand Dominant.
- H)** 91-7146-L Olson 21g Irrigating Quick Chopper, Left Hand Dominant.


**A**

**B**


**C**

### ***Irrigating Choppers***

- A)** 91-8034 91-8034 Chang 20g Irrigating Chopper, Right Curver.

### ***Lacrimal***


- B)** 91-7115 Lacrimal Cannula, 23g, Straight 25mm.  
**C)** 91-7116 Lacrimal Cannula, 23g, Straight Reinforced 20mm.  
**D)** 91-7117 Lacrimal Cannula, 23g, Curved Reinforced 20mm.  
**E)** 91-7118 Lacrimal Cannula, 19g, Side Opening, Straight, Length 28mm.




**D**

**E**

### ***LASIK, LASEK & RK***

- F)** 91-7202 Rhein LASIK Irrigation Cannula, 21g  
**G)** 91-7204 Rhein LASIK Aspiration Cannula, 19g.  
**H)** 91-8014 Slade LASIK 26g Formed Irrigation Cannula.


**F**

**G**

**H**

**A**

**B**

**C**

**LASIK, LASEK & RK**

**A)** 91-8015 Pettigrove Irrigating LASIK Cannula.

**B)** 91-8018 Banaji LASIK Irrigation Cannula, 25g.

**C)** 91-8019 Bansal LASIK Irrigation Cannula, 25g.

**D)** 91-8020 Burrato LASIK Irrigation Cannula, 25g.

**E)** 91-8021 Dishler LASIK Irrigation Cannula, 30g, 7mm.

**F)** 91-8022 Guell LASIK Irrigation Cannula, 25g.


**G)** 91-8023 Hoffman LASIK Irrigation Cannula, 25g.

**H)** 91-8024 Manche LASIK Irrigation Cannula, 25g.

**D**


**E**

**F**

**G**

**H**



**A**

**B**

**C**

### **LASIK, LASEK & RK**

**A)** 91-8025 Rubenstein LASIK Irrigation Cannula, 23g.

**B)** 91-8026 Slade LASIK Irrigation Cannula, 27g.

**C)** 91-8027 Vidaurri LASIK Irrigation Double Cannula, 25g  
91-8027-V Vidaurri LASIK Vertical Double LASIK Cannula, 25g.

**D)** 91-8028 Dishler LASIK Irrigation Cannula, 30g, 9mm




**E)** 91-8029-R Seibel LASIK Flap Irrigator & Squeegee, Right Ports


91-8029-L Seibel LASIK Flap Irrigator & Squeegee, Left Ports

**F)** 91-8030 Huang LASIK Cannula, 2 Ports, 23g

**G)** 91-8031 Plager LASIK Flap Lifter & Irrigator Cannula, Short 90 Degree Beveled Tip

**H)** 91-8033 Brown LASEK Epithelial Repositioning Cannula


**D**

**E**

**F**

**G**

**H**

**A**



**B**



**C**



**LASIK, LASEK & RK**

- A)** 91-8036 Johnston LASIK Cannula,  
Flat Tip, Rounded Distal  
Tip, 23g


**Miscellaneous**

- B)** 91-1004 Maskin Aluminum  
Handle, Male Fitting.  
**C)** 91-2000 Shepard Sideport  
Occluding Cannula.  
**D)** 91-7031 Corydon, Viscoelastic  
Irrigating Cannula, 25G,  
7/8".  
**E)** 91-7032 Viscoelastic Injection  
Cannula, 21g 7mm Tip.  
**F)** 91-7033 Kellan Viscoelastic 21g  
Cannula.


**D**




**E**



**F**






**A****B**

### **Miscellaneous**


- A)** 91-7095 Heslin Infusion Cannula,  
25g, 9mm Sandblasted Tip.  
91-7100 Heslin Infusion Cannula,  
23g, 9mm Sandblasted Tip.
- B)** 91-7102 Osher Implant Hook,  
Cannula, 24g  
With 30g Tip.
- C)** 91-7140 Gerard Anterior Chamber  
Needle Set, 23g Cannula.
- D)** 91-7205 Davis Irrigating Iris  
Retractor, 19g.


**C****D**



**A**

**B**

**C**


### *Miscellaneous*

- A)** 91-8038 Ide 0.9mm DSAEK 30g Cannula
- B)** 91-8039 Ide 1.2mm DSAEK 30g Cannula
- C)** 91-8040 Ide 1.5mm DSAEK 30g Cannula
- D)** 91-8041 Steinert DSEK Descemet Stripper w/Sandblasted Ball Tip, Irrigating
- E)** 91-8043 Friedman Irrigating IOL Manipulator
- F)** 91-8044 Phillips Lens Rotating Cannula

**D**

**E**

**F**



**A**

**B**


### ***Retrobulbar, Peribulbar & Subtenon***

- A)** 91-7078 Masket Sub-Tenon Anesthesia Cannula, 19g.
- B)** 91-7079 Rhein Open Sub-Tenon Anesthesia Cannula, 19g.
- C)** 91-7080 Atkinson Retrobulbar Cannula, 23g, 1 1/2" Long.
- 91-7085 Atkinson Retrobulbar Cannula, 25g, 1 1/2" Long.
- 91-7086 Atkinson Peribulbar Cannula, 23g, 5/8".
- D)** 91-7087 Atkinson Peribulbar Cannula, 23g, 7/8".
- 91-7088 Atkinson Peribulbar Cannula, 23g, 1 1/4".
- E)** 91-7089 Atkinson Peribulbar Cannula, 23g, 1".

**C**

**D**

**E**


<b>No.</b>	<b>Description</b>	<b>Page</b>
<b>Air Injection/Irrigation</b>		
91-1000	Delron Irrigating Handle, Two Male Fittings	2
91-1001	Delron Irrigating Handle, Male/Female Fittings	2
91-1002	Delron Irrigating Handle, Two Male Fittings, 75mm Length	2
91-1003	Aluminum Irrigating Handle, Male/Female Fitting	2
91-1100	Silicone Squeeze Bulb	2
91-1101	Bishop-Harmon Adaptor For Squeeze Bulb	2
91-1105	Bishop-Harmon Cannula, W/Bulb & Adaptor	2
91-2100	Wilbrandt-Binkhorst, 23g, Straight	2
91-2101	Wilbrandt-Binkhorst, 23g, Left	2
91-2102	Wilbrandt-Binkhorst Cannula, 23g, Right	3
91-2103	Gimbel Irrigating Cannula, 30g, 28mm	3
91-2104	Gimbel Irrigating Cannula, 27g, 25mm	3
91-2105	Gills/Welsh Cannula, 25g Angled 9mm Tip	3
91-7024	Air Injection Cannula, 25g, 5mm Tip	3
91-7025	Air Injection Cannula, 27g, 5mm Tip	3
91-7026	Air Injection Cannula, 27g, 7mm Tip	3
91-7027	Air Injection Cannula, 30g, 7mm Tip	3
91-7030	Air Injection Cannula, 30g, 5mm Tip	3
91-7035	Bishop-Harmon Cannula, 19g, 7mm Tip, 24mm	4
91-7036	Gills I/A Cannula, 25g Aspirating & 23g Irrigating	4
91-7037	Gills/Welsh Cannula, 25g, Angled 5mm Tip	4
91-7038	Gills/Welsh Cannula, 25g, Angled 8mm Tip	4
91-7039	Gills/Welsh Cannula, 25g, Angled 10mm Tip	4
91-7040	Gills/Welsh Cannula, 25g, Angled 7mm Tip	4
91-7097	Knolle Anterior Chamber Cannula, 23g, Angled 2mm Tip, 20mm	4
91-7098	Knolle Anterior Chamber Cannula, 25g, Angled 2mm Tip, 20mm	4
91-7099	Knolle Anterior Chamber Cannula, 27g, Angled 2mm Tip, 20mm	4
91-7103	Knolle Anterior Chamber Cannula, 30g, Angled 2mm Tip, 20mm	4
91-7127	Anterior Chamber Cannula, 23g Angled, 12mm	4
91-7127-S	Straight Anterior Chamber Cannula, 23g, Length 17mm	4
91-7128	Anterior Chamber Cannula, 26g Angled, 12mm	4
91-7140	Gerard Anterior Chamber Needle Set, 23g Cannula	5
91-8046	Gimbel Irrigating Cannula, Blunt Tip Angle to Side End Open 26g	21
<b>Alpha-Chymotrypsin</b>		
91-7045	Troutman Alpha-Chymotrypsin Cannula, 25g, 28mm Length	5
91-7047	Roper Alpha-Chymotrypsin Cannula, 19g	5
91-7048	Roper Alpha-Chymotrypsin Cannula, 19g, Right	5
91-7049	Roper Alpha-Chymotrypsin Cannula, 19g, Left	5
<b>Chamber Maintainers</b>		
91-7100	Heslin Infusion Cannula, 23g, 9mm Sandblasted Tip	5
91-7095	Heslin Infusion Cannula, 25g, 9mm Sandblasted Tip	5
91-8013	McIntyre III Chamber Maintainer, Cannula, 19g Tip	5
91-8037	Steinert Threaded 19g Chamber Maintainer W/Tapered Tip	5
91-8048	Gills Chamber Maintainer, 21g Tip	5

<b>No.</b>	<b>Description</b>	<b>Page</b>
------------	--------------------	-------------

**Capsule Polishers**

91-7101	Brown Capsule Polisher, Cannula, 26g	6
91-7090	Kratz Capsule Scraper, Cannula, 23g, 1mm Sandblasted Tip	6
91-7091	Knolle Capsule Scraper, 23g, 2mm Sandblasted Tip	
91-7091-1	Knolle Capsule Scraper, Cannula, 23g, 2mm Sandblasted Tip	6
91-7091-2	Lessner Capsule Scraper 27g, 3mm Tip, Sandblasted 60° Top & Bot...	6
91-7092	Jensen Capsule Scraper, Cannula, 27g, Sandblasted Tip	6
91-7093	Simcoe Capsule Scraper, Cannula, 23g, 13mm Length	6
91-7093R	Simcoe Reverse Aperture Posterior Scrubber, 23g 13mm length	6
91-7094	Drews Capsule Polisher, Cannula, 25g, Port At 6 O'Clock	7
91-7096	Kraff Capsule Polisher, Cannula, 25g Irrigation Port/Center	7
91-7101S	Brown Capsule Polisher, 26g W/Sandblasted Tip	7
91-7105	Graether Collar Button, Cannula, 23g Angled	7
91-7106	Graether Collar Button, Cannula, 23g, Straight	7
91-8002	Knolle Capsule Scraper, 30g Cannula, 2mm Sandblasted Tip	7
91-8003	Jensen Capsule Scraper, 25g Cannula, Sandblasted Tip	7
91-8004	Jensen Capsule Scraper, 27g Cannula, Semi-Sandblasted Tip	7
91-8005	Jensen Capsule Scraper, 23g Cannula, Sandblasted Tip	7
91-8006	Jensen Capsule Scraper, 25g Cannula, Semi-Sandblasted Tip	7
91-8007	Jensen Capsule Scraper 23g Cannula, Semi-Sandblasted Tip	7

**Cortex Extraction**

91-7041	Simcoe Cortex Extractor, Cannula, 23g, 0.3mm Port, 13mm Tip	8
91-7042	Simcoe Cortex Extractor, Cannula, 23g, 0.3mm Port, 21mm Tip	8
91-7043	Welsh Flat Olive Cortex Extractor, Canula, 25g, Right	8
91-7044	Welsh Flat Olive Cortex Extractor, Cannula, 25g, Left	8
91-7046	Welsh Flat Olive Cortex Extractor, Cannula, 25g, Straight	8
91-7058	Simcoe Cortex Extractor, Cannula, 23g, 0.4mm Port, 21mm	8
91-7061	Binkhorst-McIntyre, 26g, Cannula, 18mm Right	8
91-7062	Binkhorst-McIntyre, 26g, Cannula, 18mm Left	9
91-7067	Binkhorst-McIntyre, 26g, Cannula, 18mm Straight	9
91-7068	Binkhorst Cannula, 22g Right, 18mm Length	9
91-7069	Binkhorst Cannula, 22g Left, 18mm Length	9
91-7110	J-Shaped Cannula, 21g	9
91-7111	J-Shaped Cannula, 23g	9
91-7112	J-Shaped Cannula, 25g	9
91-7114	J-Shaped Cannula, 26g	9
91-8045	Anis Cortex Aspirating Cannula, Straight, 0.4mm Sideport, 23g Thin Wall	9
91-8047	Gimbel Crtex Aspirating Cannula End Open 30g	9

**Cyclodialysis**

91-7050	Kelman Cyclodialysis Cannula, 19g, Flat Right	10
91-7051	Kelman Cyclodialysis Cannula, 19g, Flat Left	10
91-7052	Randolph Cyclodialysis Cannula, 19g, Front Opening	10
91-7053	Randolph Cyclodialysis Cannula, 21g, Front Opening	10
91-7054	Gans Cyclodialysis Cannula, 19g, 32mm, Double Opening	10

<b>No.</b>	<b>Description</b>	<b>Page</b>
<b>Cystotome</b>		
91-7001	Formed Cystotome Cannula, 12mm From Tip, 22g	10
91-7002	Formed Cystotome Cannula, 12mm From Tip, 25g	10
91-7003	Formed Cystotome Cannula, 12mm From Tip, 27g	10
91-7004	Drews Cystotome Cannula, Formed, 30g, 15mm Length	10
91-7010	Kelman Cystotome Cannula, 22g Sharp, Straight	10
91-7011	Cystotome Cannula, 24g, Straight	10
91-7012	Knolle-Kelman Cystotome Cannula, 22g Long Sharp, 32mm Length	10
91-7015	Lippman Micro Cystotome Cannula, 22g, 23mm	11
91-7016	Kratz Cystotome Cannula, 22g, Angled 60 Degrees	11
91-7018	Kratz/Kelman Cystotome Cannula, 22g, Angled, 45 Degrees	11
91-7019	Kratz/Kelman Formed Cystotome Cannula, 22g	11
91-7020	Nevyas Double Sharp Cystotome Cannula, 22g, 32mm	11
91-7021	Nevyas Double Sharp Cannula, 22g, 23mm Length, Straight	11
91-7022	Kelman Cystotome Cannula, Double Sharp Front & Back, 22g	11
91-7023	Reverse Cutting Cystotome, 23g, Angled 45 Degrees	11
<b>Hydrodissection</b>		
91-7129	Hydrodissection Cannula, J-Shaped, Flattened, 30g, Straight	12
91-7131	Hydrodissection Cannula J-Shaped, Flattened, 23g, Straight	12
91-7132-A	Hydrodissection Flattened Cannula, 23g, Angled	12
91-7132-S	Hydrodissection Flattened Cannula, 23g, Straight	12
91-7133	Hydrodissection Cannula, J-Shaped, Flattened, 27g, Straight	12
91-7133-A	Hydrodissection Cannula J-Shaped 27g Flattened Angled Shaft	12
91-7134-S	Hydrodissection Cannula, 25g, Straight	12
91-7135-A	Hydrodissection Flattened Cannula, 27g, Angled	13
91-7135-S	Hydrodissection Cannula, 27g, Straight	13
91-7136-A	Hydrodissection Flattened Cannula, 30g, Angled	13
91-7136-S	Hydrodissection Flattened Cannula, 30g, Straight	13
91-7137	Hydrodissection Cannula, J-Shaped, Flattened, 25g, Straight	13
91-7138	Kellan Hydro-Elevation 27g Cannula	13
91-7141	Khoury Hydrodissection Cannula, 23g, Left	13
91-7142	Khoury Hydrodissection Cannula, 23g, Right	13
91-7206	Utrata 27g Hydrodissection "J" Cannula	13
91-8011	McIntyre III 27g Hydrodissection Cannula	14
91-8012	McIntyre III 21g Debulking Cannula	14
91-8016	Ruminson Hydrodissection Luxation Cannula	14
91-8032	Chang 27g Hydrodissection Cannula, Short 90° Beveled Tip	14
91-8032-L	Chang 27g Hydrodissection Cannula, Long 90° Beveled Tip	14
<b>I&amp;ACannulas</b>		
91-7036	Gills I/A Cannula, 25g Aspirating & 23g Irrigating	14
91-7055	Simcoe Double Barreled I/A Cannula, 23g Thin Wall	14
91-7057	Simcoe Nucleus Spatula, Cannula, 25g, Flattened Tip	14
91-7060	Simcoe Double Barreled I/A Reverse Cannula, 23g Thin Wall	14
91-7063	Gills I/A, 25g Thin Wall/Aspiration, 23g Thin Wall Irrigation	15

<b>No.</b>	<b>Description</b>	<b>Page</b>
91-7064	J-Shaped I/A Cannula, 26g Aspiration & 23g Irrigation	15
91-7071	Grandon I/A Cannula, "03", 23g X 23g, Asp. Ports Top Left	15
91-7072	Grandon I/A Cannula, "04", 23g X 23g, Asp. Ports Top Right	15
91-7122	McIntyre II 3mm Irrigating Removal Spoon, Small	15
91-7144	Khoury J-Shaped I/A Cannula, 26g Irrigation, 23g Aspiration	15
91-8042	Huang Irrigation Cannula	5
<b>Irrigating Vectis</b>		
91-7065	Knolle/Pearce Vectis Cannula, 23g, 3 Outer Ports	15
91-7066	Sheets Irrigation Vectis Cannula, 23g Without Serrations	15
91-7066-S	Sheets Irrigation Vectis Cannula, 23g With Serrations	15
91-7070	Simcoe Nucleus Loop Left, Cannula, 25g, Two Ports, 4mm	16
91-7075	Simcoe Nucleus Loop Right, Cannula, 25g, Two Ports, 4mm	16
91-7076	Lens Loop, Cannula, 3mm Wide, 3 Port	16
91-7123	McIntyre II 4.5mm Irrigating Removal Spoon, Cannula, Large	16
91-7143	Khoury 3mm Nucleus Removal Spoon, Cannula	16
<b>Irrigating Choppers</b>		
91-7145	Gillum Irrigating Phaco Splitter	16
91-7146	Olson 21g Irrigating Quick Chopper, Right Hand Dominant	16
91-7146-L	Olson 21g Irrigating Quick Chopper, Left Hand Dominant	16
91-8034	Chang 20g Irrigating Chopper, Right Curve	17
<b>Lacrimonal</b>		
91-7115	Lacrimonal Cannula, 23g, Straight 25mm	17
91-7116	Lacrimonal Cannula, 23g, Straight Reinforced 20mm	17
91-7117	Lacrimonal Cannula, 23g, Curved Reinforced 20mm	17
91-7118	Lacrimonal Cannula, 19g, Side Opening, Straight, Length 28mm	17
<b>LASIK, LASEK &amp; RK</b>		
91-7202	Rhein LASIK Irrigation Cannula, 21g	17
91-7204	Rhein LASIK Aspiration Cannula, 19g	17
91-8014	Slade LASIK 26g Formed Irrigation Cannula	17
91-8015	Pettigrove Irrigating LASIK Cannula	18
91-8018	Banaji LASIK Irrigation Cannula, 25g	18
91-8019	Bansal LASIK Irrigation Cannula, 25g	18
91-8020	Burrato LASIK Irrigation Cannula, 25g	18
91-8021	Dishler LASIK Irrigation Cannula, 30g, 7mm	18
91-8022	Guell LASIK Irrigation Cannula, 25g	18
91-8023	Hoffman LASIK Irrigation Cannula, 25g	18
91-8024	Manche LASIK Irrigation Cannula, 25g	18
91-8025	Rubenstein LASIK Irrigation Cannula, 23g	19
91-8026	Slade LASIK Irrigation Cannula, 27g	19
91-8027	Vidaurri LASIK Irrigation Double Cannula, 25g	19
91-8027-V	Vidaurri LASIK Vertical Double LASIK Cannula, 25g	19
91-8028	Dishler LASIK Irrigation Cannula, 30g, 9mm	19
91-8029-L	Seibel LASIK Flap Irrigator & Squeegee, Left Ports	19
91-8029-R	Seibel LASIK Flap Irrigator & Squeegee, Right Ports	19
91-8030	Huang LASIK Cannula, 2 Ports, 23g	19



## NOTES





**Thank You for Your  
Interest In  
Rhein Medical,  
Elevating The Quality  
Of Care In  
Ophthalmology  
Since 1988**



3360 Scherer Drive, Suite B. St.Petersburg, Florida  
800-637-4346 • Tel: 727-209-2244 • Fax: 7273418123  
Email: [Info@RheinMedical.com](mailto:Info@RheinMedical.com) • Website: [www.RheinMedical.com](http://www.RheinMedical.com)

All Rhein Products Meet The Requirements Established Within The European Union

Background Photo- Michelangelo, Pieta

A faint, stylized illustration of a person's face and upper body, rendered in a light gray tone, serving as a background for the lined paper. The illustration is positioned on the left side of the page, with the person's head tilted slightly upwards and to the right. The lines of the paper are horizontal and evenly spaced, extending across the entire width of the page.