



TYPE IIR Medical Use
SURGICAL MASK



Disposable Surgical Mask
EN14683: 2019 Type IIR
Fluid resistant ISO22609
Test Report: BST20031375202SR
Technical document: STS200313715201R
CE ECM No.1N2004051JXKH0291

100% Latex Free

Without Glass Fibres
Hypoallergenic
Easy to Breathe
Form fitting pliable nose band
High Filtration Capacity BFE ≥ 98%
Soft rounded ear-loops
Blue colour



Distributor:
www.drdoris.co.uk
info@drdoris.co.uk
GB11873332/ UAE1029

Manufacturer:
Xinxiang Kangmin Sanitary Material
Development Co., Ltd

Made in china



TYPE IIR Medical Use
SURGICAL MASK

50
pcs

MASK INSTRUCTIONS

IT'S BETTER TO USE THE RIGHT WAY TO WEAR



DISINFECTANT GEL



EAR STRAP ON
BOTH EARS



ADJUST THE BRIDGE
OF NOSE



UNFOLD THE
FOLDING PART

MATTERS NEEDING ATTENTION

1. THE MASK SHOULD BE REPLACED IN TIME. NOT RECOMMENDED FOR LONG-TERM USE
2. IT IS RECOMMENDED TO STOP USING IF THERE IS ANY MALADJUSTMENT OR AN ADVERSE REACTION DURING WEARING
3. THIS PRODUCT IS NOT WASHABLE. PLEASE MAKE SURE IT IS USED WITHIN THE VALIDITY PERIOD
4. STORE IN A DRY AND VENTILATED PLACE. AWAY FROM FIRE AND OPEN FLAMMABLES