

3 IN 1 DRUG TESTING KIT

THC 50µg/L, METH-AMPHETAMINES 300µg/L,
MDMA 500µg/L(Ecstasy)

ISO13485 ISO9001 ISO14001

ADS
AUSTRALIAN DIAGNOSTIC SERVICES



3 IN 1 DRUG TESTING KIT THC 50µg/L, METH-AMPHETAMINES 300µg/L,
MDMA 500µg/L(Ecstasy)


ADS

3 IN 1 DRUG TESTING KIT

ADS

THC 50µg/L, METH-AMPHETAMINES 300µg/L, MDMA 500µg/L(Ecstasy)

Urine drug tests screen the urine for the presence of a parent drug or its metabolites. Urine is generally considered the best medium to measure in terms of detecting drugs and drug metabolites. Compared to blood, urine offers a longer window of detection that may be one to three days for some drugs and up to 30 days for others. Urine drug testing is by far the most extensively studied and validated process for drug testing. A non-negative (positive) result indicates that the amount of drug / metabolites are above the cutoff limits for the test.

Substance	Cut off µg/ml	Detectable in Urine	Detection Period (days)	Contents
MethAmphetamine (MET)	300 µg/ml	2 - 5 hours after use	2 - 4	1x Multi Panel Urine Drug Testing Kit 1x Instruction manual
Marijuana (THC)	50 µg/ml	2 - 5 hours after use	15 - 30	ISO13485 ISO9001 ISO14001 Approved
Ecstasy (MDMA)	500 µg/ml	2 - 5 hours after use	2 - 3	 0 123456 789012

info@diagnosticservicesaustralia.com

ADS

3 IN 1 DRUG TESTING KIT

ADS

3 IN 1 DRUG TESTING KIT THC 50µg/L, METH-AMPHETAMINES 300µg/L,
MDMA 500µg/L(Ecstasy)