

# Aqcura™ SHOWER CHAIR

## PRODUCT SPECIFICATION SHEET

27<sup>th</sup> May 2020

**Name:** Aqcura™ Shower Chair or Shower Commode.



**Model Nos:** SC460PBR with Padded Back Rest & SC460SBR with 3-Slat Back Rest. **Model Nos:** SC460/O with Open front seat & SC460/C with Closed front seat.

**TGA Number:** 335975 on the Australian Register of Therapeutic Goods (ARTG.)

**Testing Background:** The Aqcura Shower Chair has been tested to Australian and New Zealand Standards for Shc with 1000 hour exposure in salt spray solution. AS/NZS 3973:2009 Type B

The Aqcura™ Shower Chair was tested and received a Certificate of Compliance with NATA Report No. 20-2902968. Complies with the requirements of AS/NZS 3973:2009 for maximum user weight at 200kg.

**Product Features:** Tubular 304 grade stainless steel construction with moulded seat base, back rest and arm rests. This is an attendant propellant chair with sliding weight bearing foot rest. Equipped with four swivel locking castors. Castor housing made from high grade synthetic materials with stainless steel metal parts.



Dimensional drawing of Aqcura™ Shower Chair

SIDE VIEW

FRONT VIEW

ACCURA. SHOWER CHAIR:  
OVER TOILET FITTING PROBLEM

SC460 SHOWER CHAIR

